

Comparison of Various models of Stand

SN	Particulars			
MODEL		SS-01	SS-02	SS-03
		(NORMAL DUTY)	(MEDIUM DUTY)	(HEAVY DUTY)
1	Shape of the Base	C	C	C
2	Weight of the Cast Iron Base	6 Kg.	13 Kg.	19 Kg.
3	Diameter of SS-316 Stand Rod	16 mm	25 mm	25 mm
4	Length of the Stand Rod	1 Meter	1 Meter	1 Meter
5	Duty	N-(Normal)	M-(Medium)	H-(Heavy)